


# Adria Multi-Purpose Drainage Catheter



The material of Adria's Multi-Purpose drainage catheter enhances patient comfort with a smooth, Kink-resistant flexible shaft.

Large drainage holes and generous lumen allow for optimal flow dynamics.

Slick surface and smooth, tapered tip facilitates catheter insertion and final placement.

The Adria Multi-Purpose Drainage Catheter is an excellent choice for percutaneous nephrostomy, abscess drainage, and general fluid collections.

The Adria Multi-Purpose Drainage Catheter can be introduced according to either the "one step", or Seldinger Technique (coaxial or parallel access puncture).

Available in J-Tip, Pigtail, and locking pigtail configuration.

## Catalogue Codes

80106	6 Fr. x 30 cm, J-Tip
80306	6 Fr. x 30 cm, Pigtail
80506	6 Fr. x 30 cm, Locking Pigtail
80107	7 Fr. x 30 cm, J-Tip
80307	7 Fr. x 30 cm, Pigtail
80507	7 Fr. x 30 cm, Locking Pigtail
80108	8 Fr. x 30 cm, J-Tip
80308	8 Fr. x 30 cm, Pigtail
80508	8 Fr. x 30 cm, Locking Pigtail
80109	9 Fr. x 30 cm, J-Tip
80309	9 Fr. x 30 cm, Pigtail
80509	9 Fr. x 30 cm, Locking Pigtail
80110	10 Fr. x 30 cm, J-Tip
80310	10 Fr. x 30 cm, Pigtail
80510	10 Fr. x 30 cm, Locking Pigtail
80112	12 Fr. x 30 cm, J-Tip
80312	12 Fr. x 30 cm, Pigtail
80512	12 Fr. x 30 cm, Locking Pigtail

Packaged 5 per carton