

Hysterosalpingography Catheter


This catheter allows an easy and atraumatic injection of saline solution into the uterus during trans-vaginal echography.

The special shape of the balloon has been designed to avoid saline solution's reflux and to perform an appropriate dilation.

The balloon is produced in Latex to prevent the risk of allergies.

Main Features

- ✓ Radiopaque catheter
- ✓ Flexible and transparent balloon
- ✓ Round tip with side holes to guarantee an outstanding contrast disposal



Catalogue Codes

Order Code	Description
76005	05Fr x 30cm
76007	07Fr x 30cm

Packaged 5 per carton

ADRIA s.r.l.
40033 Casalecchio di Reno (BO)
Via della Bastia, 9 - ITALY
Tel. 051.6199920
Fax 051.6183098
E-mail: info@adriamedical.com
www.adriamedical.com

ADRIA s.r.l.