



Adria "Illinois Type" Sternal/Iliac Aspiration Needle



Explant Needle

This single-use, adjustable-length needle is specifically designed to obtain bone marrow samples by aspiration.

Penetrates easily into the bone, ensuring safe and simple aspiration of marrow from the iliac crest or sternum.

The adjustable depth guard controls the depth of entry into the bone.

The stainless steel luer lock ensures secure syringe attachment.

Catalogue Codes

64115	15 ga x 10-44 mm
64315	15 ga x 25-60 mm
64116	16 ga x 10-44 mm
64316	16 ga x 25-60 mm
64118	18 ga x 10-44 mm

Packaged 10 per carton

Explant

64114	14 ga x 37-71 mm
-------	------------------

Packaged 5 per carton