

Iso-Sure Breast Localization Needle

Localization needles for non-clinically palpable mammary lesions.

The hub design allows a quick and smooth entry of the needle, while the hook provides optimal anchoring strength.

- Echogenic Tip
- Graduated centimeter markings
- Movable stopper for precise placement
- Standard 20 gauge needle

Catalogue Codes

Iso-Sure, Luer lock grip

41315	3 cm needle, 13 cm wire
41515	5 cm needle, 15 cm wire
41720	7 cm needle, 17 cm wire
41920	9 cm needle, 19 cm wire
41435	15 cm needle, 25 cm wire

Packaged 10 per carton

Iso-Sure, Fine grip

42315	3 cm needle, 13 cm wire
42515	5 cm needle, 15 cm wire
42720	7 cm needle, 17 cm wire
42920	9 cm needle, 19 cm wire
42435	15 cm needle, 25 cm wire

Packaged 10 per carton

