


Nephrostomy High Pressure Balloon Dilatation Catheter



This catheter is designed for the dilatation of nephrostomy tract.

It guarantees an easy placement of the working sheath.

Main Features:

- ✓ 7 Fr catheter with high pressure stopcock
- ✓ Two radiopaque markers
- ✓ PTFE radiopaque 17 cm length sheath
- ✓ 12 cm balloon length
- ✓ 10 mm inflated balloon with uniform pressure distribution (30 Fr)
- ✓ Maximum operating pressure: 17 ATM

Order Code	Catheter Size	Balloon
25312	7Fr x 55cm	10mm x 12cm

Packaged 5 per carton