


# Paracentesis kit



This single-use, retractable tip needle is specifically designed to obtain safe drainage thanks to the possibility to avoid the most sensitive tissue.

Echogenic Tip enhances visualization during ultrasound procedures.

Graduated centimeter markings ensure precise placement.

The Kit includes:

- 8000ml drainage bag with non return valve, drain cock and hanging.
- 60ml plastic syringe luer lock.
- radiopaque trocar point cannula with side hole and retractable stylet.

## Catalogue Codes

72015	15 ga x 10 cm
72018	18 ga x 5 cm
72020	18 ga x 10 cm

*Packaged 10 per carton*