

Angioplasty, PTA Balloon Catheter Adapted to .018i and .035i



OUTSTANDING TRACKABILITY

Short tapered tip design

Highly kink-resistant shaft design

Low deflated profile;

Excellent pushability

EXCELLENT DILATION RESULTS

Uniform pressure distribution;

Highly scratch and puncture resistant

High rated burst pressure

How to order

Angioplasty PTA Balloon Catheter .035i

Catalogue N° Usable Length (75 cm)	Catalogue N° Usable Length (120 cm)	Balloon Diameter (mm)	Balloon Length (mm)	Introducer (F)	Rated Burst Pressure (atm)	Recommended Guide wire (in)
77004-40200	77005-40200	4	20	5	20	.035
77004-40400	77005-40400	4	40	5	20	.035
77004-40600	77005-40600	4	60	5	20	.035
77004-41000	77005-41000	4	100	5	18	.035
77004-50200	77005-50200	5	20	5	20	.035
77004-50400	77005-50400	5	40	5	20	.035
77004-50600	77005-50600	5	60	5	18	.035
77004-50800	77005-50800	5	80	5	18	.035
77004-60200	77005-60200	6	20	5	19	.035
77004-60400	77005-60400	6	40	5	19	.035
77004-60600	77005-60600	6	60	5	18	.035
77004-60800	77005-60800	6	80	6	17	.035
77004-70200	77005-70200	7	20	5	16	.035
77004-70400	77005-70400	7	40	5	16	.035
77004-70600	77005-70600	7	60	6	17	.035
77004-70800	77005-70800	7	80	6	17	.035
77004-80200	77005-80200	8	20	5	15	.035
77004-80400	77005-80400	8	40	5	15	.035
77004-80600	77005-80600	8	60	6	15	.035
77004-90200	77005-90200	9	20	6	15	.035
77004-90400	77005-90400	9	40	6	15	.035
77004-90600	77005-90600	9	60	6	15	.035
77004-10400	77005-10400	10	40	6	13	.035
77004-10600	77005-10600	10	60	7	14	.035
77004-12200	-	12	20	7	12	.035
77004-12400	-	12	40	7	12	.035
77004-14200	-	14	20	7	10	.035
77004-14400	-	14	40	7	10	.035

Scatola/Box: 1 Unità/unit

Angioplasty PTA Balloon catheter .018i

Catalogue N° Usable Length (90 cm)	Catalogue N° Usable Length (120 cm)	Balloon Diameter (mm)	Balloon Length (mm)	Introducer (F)	Rated Burst Pressure (atm)	Recommended Guide wire (inch)
-	77007-20200	2,0	20	4	20	.018
-	77007-20400	2,0	40	4	20	.018
-	77007-25200	2,5	20	4	20	.018
-	77007-25400	2,5	40	4	20	.018
-	77007-30200	3,0	20	4	20	.018
-	77007-30400	3,0	40	4	20	.018
-	77007-31000	3,0	100	4	18	.018
77006-40200	77007-40200	4,0	20	4	20	.018
77006-40400	-	4,0	40	4	20	.018
77006-40600	-	4,0	60	4	18	.018
77006-50200	-	5,0	20	4	20	.018
77006-50400	-	5,0	40	4	20	.035
77006-50800	-	5,0	80	4	16	.018
77006-60200	-	6,0	20	4	17	.018
77006-60400	-	6,0	40	4	17	.018
77006-60600	-	6,0	60	4	14	.018
77006-70200	-	7,0	20	4	14	.018
77006-70400	-	7,0	40	4	14	.018
77006-70600	-	7,0	60	4	14	.018
77006-80200	-	8,0	20	5	14	.018
77006-80400	-	8,0	40	5	14	.018
77006-90200	-	9,0	20	5	13	.018
77006-90400	-	9,0	40	5	13	.018

Scatola/Box: 1 Unità/unit