


Retrieve® Disposable Bone Marrow Aspiration Needle



The Retrieve®, bone marrow biopsy needle, was specifically designed for the efficient and effective sampling of bone marrow from the iliac crest.

Razor sharp cutting edges, coupled with sturdy construction and design, allows the physician to sample with reduced effort and heightened confidence.

Tapered cannula ensures tissue retention while leaving the cellular structure intact.

Single patient use reduces costly and time consuming clean up and eliminates the risk of cross contamination.

Available in both standard bone marrow biopsy and transplant versions.

Catalogue Codes

Standard

50001	7 ga x 10 cm
50002	7 ga x 15 cm
50003	8 ga x 10 cm
50004	8 ga x 15 cm
50005	9 ga x 10 cm
50006	9 ga x 15 cm
50007	11 ga x 10 cm
50008	11 ga x 15 cm
50009	13 ga x 10 cm

Packaged 5 per carton

Transplant

51001	11 ga x 9 cm
51002	11 ga x 13 cm
51003	13 ga x 7 cm
51004	13 ga x 11 cm
51005	15 ga x 7 cm
51006	15 ga x 11 cm

Packaged 5 per carton