


Adriastone Stone retrieval Catheter



Adristone is specially designed for treatment of ureter stones.

Its excellent flexibility in conjunction with its kink resistant material guarantees an atraumatic access.

The retrieval of the stones results very easy thanks to helical and round shape of the wires.

The package includes:

- ✓ Handgrip and syringe port for one hand procedure
- ✓ Basket
- ✓ Shaft

Radiopaque

Latex-free

Order Code	Description
23202	2.5Fr x 90cm x 4W
23302	3.0Fr x 90cm x 4W
23402	3.5Fr x 90cm x 4W
23502	4.0Fr x 90cm x 4W
23602	5.5Fr x 90cm x 4W
23702	5.0Fr x 90cm x 4W

Packaged 5 per carton