



# Vertebroplasty Set



## The Adria Vertebroplasty Set includes:

- ✓ An Injector Gun
- ✓ A special Vertebroplasty Needle
- ✓ One unit of radiopaque cement
- ✓ A connection tube between gun and needle



## Vertebroplasty Injector Gun:

- ✓ The large screw allows a rapid injection assuring a rapid and precise performance
- ✓ The ergonomic design guarantees a comfortable utilization during the procedure
- ✓ The safety button permit an instant pressure relief in case of leakage, stopping the cement injection

## Special Vertebroplasty Needle:

- ✓ The small, non-metal handgrip reduces artefacts under CT
- ✓ The small handgrip's wings facilitate the insertion and the removal of the needle, allowing its rotation into the vertebral body
- ✓ Easy insertion is also guaranteed by the special tip sharpening
- ✓ The bayonet luer lock ensures the stability of the stylet into the cannula during the procedure

## Catalogue Codes

### VERTEBROPLASTY SET

Order Code	Description
49009120	09ga x 12cm
49009150	09ga x 15cm
49011120	11ga x 12cm
49011150	11ga x 15cm
49013120	13ga x 12cm
49013150	13ga x 15cm

Packaged 1 per carton

### VERTEBROPLASTY NEEDLE

Order Code	Description
49112	09ga x 12cm
49115	09ga x 15cm
49212	11ga x 12cm
49215	11ga x 15cm
49312	13ga x 12cm
49315	13ga x 15cm

Packaged 5 per carton

### VERTEBROPLASTY CEMENT

Order Code	Volume
49100020	20ml

Packaged 1 per carton

### VERTEBROPLASTY HAMMER

Order Code	Weight
49218524	240g

Packaged 1 per carton